



*Magic on your face ...*

### CLIKINS Orenu

It marks the color pigments inhale into the needle by injecting them to skin by a back-and-forth motion MTS machine helps medical supplies be absorbed to skin by irritating.



Manufacturer: Biz&Mold Co.,Ltd  
ADD : B-401,Songdo Smartvalley,30,Songdomirae-ro,  
Yeonsu-gu,Incheon,Korea  
T: +82-32-876-0919 F: +82-32-232-0924 M:+82-10-2367-7997  
E: biznmold@hanmail.net H : http://clikins.com

Magic on your face...

<http://clikins.com>



### CLIKINS JUST 6 DAYS

Cosmetic absorbing devices that help improve cosmetics and maximize the effects of cosmetics on cosmetics

CLIKINS roller is an innovative skin roller system consisting of micro-needle, LED, and vibration massage. The small spits on this small tool create numerous micro spots on the skin to activate the skin penetration of the skin care product. CLIKINS rollers use a very special intelligent chip to control the wavelength of light and ultra-low vibration. Clikins Roller ha helps to absorb the absorption ratio of 60 to 70 % (20 to 30%).



#### Home care set components

The best avengers in the category of cost-effectiveness in one home care set!


It is the best home care program that can start with the effective skin absorption of the pig-placenta ample, and can easily proceed to the skin calming mask pack at once.

### CIESKIN

The CIESKIN program is the starting program for beautiful skin in Korea. It is the best program to apply to skin wrinkles, whitening, lifting and skin signs.

Beautiful skin miracle only for you!

Using micro current and 650nm wavelength, beauty device helps to absorb the ample.



#### CLESKIN PROGRAM

- 1 Clean your skin (cleansing)
- 2 Relax with warm-water
- 3 Soft skin massage
- 4 Applying solutions
- 5 Apply the solution from the face to the skin with a soft stamp
- 6 Apply all remaining solution to face and put mask pack
- 7 Mask the pack and once more MTS Toning and Cool Soothing
- 8 Finish after application of regenerated cream and sunscreen

### CLIKINS H + COLOR PIGMENT

#### 1. Object

Apparatus for applying a Micro-needle Therapy by vibration.

#### 2. Performance

A device for puncturing skin-penetrating drugs or cosmetics supports with a micro-needle guided between extended and retracted positions.

#### 3. General maintenance

- Keep it with its power off. Wipe its surface out with a dry cloth.
- Keep it away from direct sunlight or high temperature and humidity
- Do not leave the measurement, repair or alter of the device to non-professional staff.